

ПРЕДВАРИТЕЛЬНЫЕ ФИЛЬТРЫ MILLIGARD



Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97

Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

<https://millipore.nt-rt.ru> || <mailto:mer@nt-rt.ru>

Milligard® Filters

Superior filters for critical prefiltration applications

Milligard filters provide the finest prefiltration available. Milligard filters are constructed of homogeneous polymers of cellulose esters on a polyester web resulting in a versatile filter with high reproducibility and efficacy due to the highly retentive, non-fiber releasing membrane. The homogenous cellulose ester membrane provides the high retention needed during critical prefiltration steps while protecting more expensive downstream filtration devices from premature plugging.

Milligard filters are available in single and double layer formats.

Benefits

- Effectively protect expensive downstream filters from premature plugging
- High dirt holding capacity combined with high retention efficiency
- Unparalleled thermal and hydraulic stress resistance
- Validated to withstand multiple sterilization cycles
- Ideal for designing scalable solutions from bench top to full-scale manufacturing



Media Types

Milligard
0.2 µm
0.5 µm
1.2 µm
0.5/0.2 µm
1.2/0.5 µm
Milligard pore sizes are nominal

Filter Formats

- OptiScale® small scale disposable capsule filters
- Opticap® XL 2 disposable capsule filters

APPLICATIONS

Cell Culture Media

Effectively remove particles and colloidal contaminants without obstructing the flow of vital media constituents. Milligard prefilters provide the high flow rates and consistently high throughput required to accomplish the separation in a timely and cost-effective manner.

Final Aseptic Fill of Small Volume Parenterals (SVP)

Extend the service life of downstream sterilizing filters by removing colloidal and particulate contaminants.

Large Volume Parenterals (LVP)

Extend the service life of downstream sterilizing filters by removing colloidal and particulate contaminants. Extremely robust and able to withstand high operating pressures, Milligard filters offer dependable retention and high flow rates needed for high volume manufacturing.

Plasma Proteins—Human Albumin

Remove colloids, aggregated and non-product proteins, lipids and particles prior to downstream purification without retaining the fractions of interest. Milligard filters provide superior protection of downstream sterile filters and/or chromatography columns.

Serum

Effectively remove lipids, colloids, and particles from serum before the final sterilizing filtration without obstructing the passage of serum proteins. Milligard filters provide high flow rates and consistently high throughput.

REGULATORY COMPLIANCE

Milligard filters are designed, developed, and manufactured in accordance with a Quality Management System approved by an accredited registering body to an ISO® 9000 Quality Systems Standard and are shipped with a Certificate of Quality. Each Opticap XL 2 capsule filter is supported by a Validation Guide for compliance with regulatory requirements.

For traceability and easy identification, each filter is marked with identifying characteristics.

MULTIPLE FORMATS AVAILABLE

Milligard filters are available in single or double layer, three formats, many pore sizes, and multiple configurations that vary by filtration area and the type of inlet/outlet connection.

OPTISCALE PROCESS DEVELOPMENT SCREENING TOOL

OptiScale disposable capsule filters with Milligard media provide a convenient small-volume option for process screening and scaling. These "drop in" filters are ideal for evaluating biopharmaceuticals. OptiScale capsule filters offer speed-to-market strategies for efficiently developing compounds and biotherapeutics.

The OptiScale disposable capsule is ideally suited for process development and screening. OptiScale capsules are faster and easier to set up than conventional 47 mm discs.



OPTICAP XL 2 DISPOSABLE CAPSULE FILTERS

Convenient and Easy to Use

Opticap XL 2 capsule filter eliminates the time and expense associated with assembling, cleaning, and validating stainless steel housings. Adjustable, easy-to-turn, upstream vents and drain valves with O-ring seals and hose barb connections allow for easy process control. Other ease-of-use features include flow direction arrows and ribbed housing for easy gripping even with gloved hands.

The Right Size

The Opticap capsule product family provides a wide range of filtration areas to fit all of your application needs and to allow easy scale-up of your small volume filtration steps to larger, full-scale filtration processes.

The Right Connections

Self-contained and disposable, Opticap XL 2 capsule filters are supplied with a choice of inlet and outlet connections to optimize your filtration process, including sanitary flanges which provide a high flow rate, fractional sanitary flanges, and hose barbs.

Robust Construction

Opticap XL 2 capsule's design allows unparalleled thermal and hydraulic stress resistance in a disposable filter, resulting in reliability, high confidence in the sterility process, and improved cleanliness.

SPECIFICATIONS (OptiScale Capsule Filters)

Nominal Dimensions	
Maximum length:	82 mm (3.24 in.) with flange inlet/hose barb outlet 74 mm (2.91 in.) with flange inlet/flange outlet 94 mm (3.70 in.) with hose barb inlet/hose barb outlet
Diameter:	69 mm (2.75 in.)
Weight:	2.3 oz (67 g)
Filtration Area	13.8 cm ²
Materials of Construction	
Filter media:	Mixed esters of cellulose
Structural components:	Polycarbonate
Vent cap:	PVDF
Internal seal rings:	Viton® fluoroelastomers
Housing Vent	Adjustable vent with male luer and female Luer-Lok™ connections on inlet side of device.
Maximum Inlet Pressure	5.5 bar (80 psi) at 25 °C
Oxidizable Substances	Capsules meet the requirements of the USP Oxidizable Substance for Sterile Water for Filtration Test after a water flush of ≤100 mL.
Sterilization	May be autoclaved for 3 cycles of 60 minutes at 121 °C.
Good Manufacturing Practices	These products are manufactured in a Millipore facility which adheres to FDA Good
Practices	Manufacturing Practices.
Component Material	Component materials meet the requirements of the USP <88> Reactivity Tests for Class VI Plastics.
Toxicity	This product is non-toxic per the USP <88> Safety Test.
Indirect Food Additive	All components meet the FDA Indirect Food Additive requirements cited in 21 CFR 177–182.


SPECIFICATIONS (Opticap XL 2 Capsule Filters)

Nominal Dimensions	
Maximum length:	14.2 cm (5.6 in.)
Nominal Diameters	
Body diameter:	8.4 cm (3.3 in.)
Vent to vent diameter:	12.5 cm (4.9 in.)
Filtration Area	
Single layer:	0.10 m ² (1.1 ft ²)
Double layer:	0.07 m ² (0.8 ft ²)
Materials of Construction	
Milligard Media:	Mixed esters of cellulose
Supports:	Polyester
Structural components*:	Polypropylene
Vent O-rings:	Silicone
Vent/Drain	¼ in. hose barb with double O-ring seal
Maximum Inlet Pressure	
	5.5 bar (80 psi) at 23 °C 2.8 bar (40 psi) at 60 °C 1.0 bar (15 psi) at 80 °C
Maximum Differential Pressure	
Forward:	3.4 bar (50 psid) at ambient room temperature.
NVR Gravimetric Extractables	
	After autoclaving and a 24 hour soak in ASTM® Type 1 reagent grade water at controlled room temperature:
Milligard Media:	≤ 15 mg
After a flush of:	1000 mL
Oxidizable Substances	
	Meets the requirements of the USP Oxidizable Substances Test after a water flush of ≤ 1000 mL
Bacterial Endotoxin	
	Aqueous extraction contains < 0.5 EU/mL as determined by the Limulus Amebocyte Lysate (LAL) Test.
Sterilization	
	May be autoclaved for 3 cycles of 30 minutes at 121 °C. (Cannot be steam sterilized in-line.)
Good Manufacturing Practices	
	These products are manufactured in a Millipore facility which adheres to FDA Good Manufacturing Practices.
Non-Fiber Releasing	
	Milligard media meets the criteria for a "non-fiber releasing" filter as defined in 21 CFR 210.3 (b) (6).
Component Material Toxicity	
	Component materials were tested and meet the criteria of the USP <88> Reactivity Test for Class VI Plastics. Milligard filters meet the requirements of the USP <88> Safety Test utilizing a 0.9% sodium chloride extraction.
Indirect Food Additive	
	All component materials meet the FDA Indirect Food Additive requirements cited in 21 CFR 177–182.


*Cage, core, end caps and capsule housing

TYPICAL CLEAN WATER FLOW RATES


Opticap XL 2 Capsule with Milligard
Media 0.2 μm Nominal (KW03)




Opticap XL 2 Capsule with Milligard
Media 0.5 μm Nominal (KW06)




Opticap XL 2 Capsule with Milligard
Media 1.2 μm Nominal (KW19)



Opticap XL 2 Capsule with Milligard
Media 0.5/0.2 μm Nominal (KWSS)



Opticap XL 2 Capsule with Milligard
Media 1.2/0.5 μm Nominal (KWSC)



Opticap XL Capsule Legends

Refer to Connection Type


TT = 38 mm (1 1/2 in.) Sanitary Flange Inlet and Outlet

FF = 19 mm (3/4 in.) Sanitary Flange Inlet and Outlet


HH = 14 mm (1/2 in.) Hose Barb Inlet and Outlet

ORDERING INFORMATION


OptiScale Capsule Filters



Opticap XL 2 Capsule Filters



Cartridge Filters





Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97

Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

<https://millipore.nt-rt.ru> || <mailto:mer@nt-rt.ru>