

ХРОМАТОГРАФИЧЕСКИЕ СИСТЕМЫ MOBIUS CHROM 20



Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97

Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

<https://millipore.nt-rt.ru> || <mailto:mer@nt-rt.ru>

Mobius® Chrom 20 System

For Purification of Clinical and Process Scale Biologics

The Mobius® Chrom 20 system is a flexible, automated, single-use chromatography system that enables consistent and reliable separation and purification of mAbs, vaccines, plasma, and therapeutic proteins at clinical and process scales.

The innovative single-use flowpath reduces contamination risk and gives you the confidence to tackle any chromatography challenge, while the flexible design featuring a standard hardware platform that can be combined with a selection of options, empowers you to meet your unique process requirements.

The Mobius® Chrom 20 system delivers optimal operational flexibility for processes up to 2000L with flow rates up to 20L/min, a 5:1 turn down ratio at process scale, and enhanced isocratic and linear gradient accuracy.



Benefits

Easy to Use

Our unique clamshell and single-use Flexware® clamshell assemblies were designed to offer ultimate ease of use from preparation to production. Simply insert the assembly into the clamshell and connect the tubings.

Safe & Reliable Process

Our 100% single-use, gamma-irradiated Flexware® assemblies provide quick product changeover with reduced contamination risk and cleaning for greater efficiency and process security.

The high gradient accuracy of the system enhances separation performance making it ideal for challenging processes.

Flexible

Easily scale-down to lower column sizes while maintaining high accuracy to meet broader processing requirements. Ability to cover multiple scales and applications offers maximum return on investment.

User-friendly

The system was designed with operator comfort and safety in mind and features a movable screen for operator convenience. The stainless-steel surface is easy to clean and the accessible components facilitate maintenance.

System Components

The Mobius® Chrom 20 system is made up of a base system and configurable options that can be added to design a system that best meets your purification process and challenges.

The base cart is the heart of the system and houses the clamshell and automation software.


The unique clamshell design ensures simple installation of the flowpaths so your process is up and running rapidly, while the optional clamshell lift enables a single technician to assemble and maintain the system.

The system also shares common elements with our other unit operations, such as an intuitive software interface, assembly material (PureFlex™ film), and maintenance components, making training and validation easier across scales and applications.

Mobius® Chrom 20 System Configurations

The basic configuration includes:

- Clamshell
- Pump cart
- Base cart
- Post-column instrumentation



The system options and accessories include:



Pre-column 0,2µm filter and its support for protecting the column



Pre-column instrumentation for conductivity and pH measurement



Keyboard support for enhanced ease of use



Clamshell lift for one-man assembly and maintenance

Mobius® Chrom 20 Flexware® Assemblies

Single-use Flexware® assemblies deliver equivalent performance to that of traditional stainless-steel systems, while providing maximized contamination control. Assemblies minimize dead space, decreasing the risk of contamination to virtually zero, while ensuring high gradient accuracy. The assemblies are delivered gamma-irradiated, ready-to-use, and are clearly labeled, significantly reducing the risk of installation error and the challenges typically associated with installing single-use flow paths.

Benefits of our 100% Single-use Flowpaths

- Fast changeover and reduced contamination risk from previous batches in multi-product and multi-scale productions
- Combination of clamshell and Flexware® clamshell assembly offers easy installation compared to classic tubing
- Reduced hold-up volume contributes to the enhanced accuracy of the system
- Reduced number of connections versus traditional single-use assemblies decreases the risk of integrity failure
- The clamshell assemblies are 100% integrity tested during manufacturing, ensuring reliable performance



Flexware® clamshell assembly labels



Mobius® Chrom 20 system clamshell and Flexware® clamshell assembly

Optimized, Robust Chromatography Performance

Exceptional gradient accuracy performance

The system is designed to deliver high accuracy gradients using pump-based mixing over an extended flowrate range. The test results below illustrate the better than 2% accuracy obtained from 4 to 20 L/min over a mixing range of 10-90% in step and linear gradient modes. This unique performance enables enhanced separation in a fully single-use mode to meet any chromatography challenge.




Figure 1: Linear Gradient (%) at 4 L/min, 12 L/min and 20 L/min




Figure 2: Step Gradient (%) at 4 L/min, 12 L/min and 20 L/min

Low pressure drop over full flow rate range

The system is designed to minimize pressure drop across the full flow rate range. The system pressure drop is below 1 bar at max flowrate.




Figure 3: Pressure drop (bar) at 4 L/min, 12 L/min and 20 L/min

Mobius® Chrom 20 System and Flexware® Assemblies Specifications

General System Specifications

Specification		Mobius® Chrom 20 System
Flow Range		In gradient mode within accuracy specifications: 240-1200L/h In Isocratic mode: 24-1200L/h
Pumps	Model	Quattroflow 1200S x 2
	Maximum Pressure	4 bar
Valves	On/Off Valves	Inside Clamshell: Automatic pneumatic piston
		Inlet Manifold: Automatic pinch valve for 5/8" ID tubing
	Valve control (after Bubble Trap)	Pneumatic piston with pressure regulator

Mechanical Specifications: System and Assemblies

Specification		Mobius® Chrom 20 System
System Dimensions	Pump & Base carts connected H x W x D in mm	1727 x 1715 x 800 ± 20 mm (68 x 67.5 x 31.5 in.)
Net Weight	Base cart with Clamshell	365 kg (805 lbs)
	Pump cart with 2 Pumps and Manifold	160 kg (353 lbs)
Materials of Construction	Wetted Components USP<88> Class VI compliant material	
	Tubing	Silicone
	Flexware® Clamshell Assembly	Pureflex™ Film, polyethylene fittings and silicone gaskets, polypropylene
	Pump Head	EPDM, Santoprene® and polypropylene
	Connectors	Polypropylene, polysulfone
	Flowmeter	Polysulfone and Hastelloy C22 electrodes
	UV, Conductivity and pH Sensor (MU)	Quartz, EPDM and stainless steel 316L
	UV, Conductivity and pH Sensor (SUC)	Polysulfone Quartz, EPDM and stainless steel 316L (pins only)
	Non-Wetted Components	
	Carts frame & panels	Stainless steel 304L
	Clamshell (rear)	POM
	Clamshell (front)	PMMA
	Valve Pads	Silicone
Casters		8 (4 free, 4 locking) Inner : Polyamide 6 Rolling area : Polyurethane
Connections	Inlets (10)	¾" TC
	Outlets, 3 Fractions + 1 waste	¾" TC
	Filters	¾" TC
	Column	¾" TC
	Pneumatic	¾" TC

Operating Specifications

Specification	Mobius® Chrom 20 System	
Total Operating Time	Flexware® Assemblies	Validated for 24 hours operation 50 cycles (50 valve openings and 50 valve closings) per valve The integrity of the installation may be compromised if used for more than 50 cycles
	Valve Pads	Do not exceed 1200 cycles or 6 months of operation
Product Temperature Range		20 to 30°C (68 to 86°F)
Maximum Pressure	Pump 1 Manifold	0 to 2 bar (0 to 29 psi)
	Pump 2 Manifold	0 to 2 bar (0 to 29 psi)
	Pump Assembly	0 to 4 bar (0 to 58 psi)
	Bubble Trap Assembly	0 to 4 bar (0 to 58 psi)
	Clamshell assembly	0 to 4 bar (0 to 58 psi)
	Precolumn Filter Assembly	0 to 4 bar (0 to 58 psi)
	Pre-column Instrumentation Assembly	0 to 4 bar (0 to 58 psi)
	Post-column Instrumentation Assembly	0 to 2 bar (0 to 29 psi)
	Column Assembly	0 to 4 bar (0 to 58 psi)
System Operating Temperature		20 to 30°C (68 to 86°F)
Operating Humidity		10 to 90% (non condensing)

Power and Air Supply Specifications

Specification	Mobius® Chrom 20 System
Power Supply	220-240VAC, 50/60 Hz, 1 phase, 6.5 A max or 100-120VAC, 50/60 Hz, 1 phase, 13.6 A max
	Maximum consumption 1.6 KW
Pneumatic Air Supply	6 bar (87 psi) min, 10 bar (145 psi) maximum, oil free
	Maximum consumption 4 L/min

Instrument Specifications and Performances

Specification	Tag*	Range/Setting/Type/Accuracy	Process Connection
Pressure Indicators	PIT001 PIT002 PT003	0-4 bar (0-58 psi) ± 0.2 bar (2.9 psi)	Non-intrusive
High Pressure Switches	PSH001 PSH002	4.2 bar (60 psi)	Non-intrusive
Single-Use Magnetic Flowmeters	FT001 FT002	± 2% MV between 0,4 L/min and 20 L/min	In-line
Temperature Sensors	TE201	2-30°C (35.6°F to 86°F) ± 1°C (33.8°F)	Combined with conductivity sensor
Single-use Conductivity Sensors	AE102 AE201A AE201B	Process: 0-100 mS/cm ±2% FS Cleaning: 0-50 µS/cm ±2% FS	In-line
pH Sensors	AE103 AE202	1-14 pH ±0.1 pH	In-line
Air Sensors	XS001 XS002	N/A	Non-intrusive
UV Sensors	AE203/204	AF45 : 0-3 AU ±2% FS OPL: 1 or 2.5 mm Wavelength: 280 nm AF46 : 0-2 AU ±2% FS OPL: 1 or 2.5 mm Wavelength: Dual 254/280 nm or 280/300 nm	In-line

Do not use UV Sensors in a condensing atmosphere. Condensation may lead to erroneous sensor readings.

*Refers to the labels on the system hardware

Ordering Information

Hardware

Catalogue Number	Description	Base Cart	Bubble Trap Support	Pump Cart with 2 Pumps and 2 Single-use Flowmeters	Single-use Pre-column Instrumentation Holder	280 nm	280/300 nm	254/280 nm
		Standard	SU Post-column Instrumentation					
China								
MCPXM20CNSB01	Mobius® Chrom 20 Post SU 280 nm	Included				●		
MCPXM20CNSB02	Mobius® Chrom 20 Post SU 280-300nm						●	
MCPXM20CNSB03	Mobius® Chrom 20 Post SU 254-280nm							●
MCPXM20CNSBS1	Mobius® Chrom 20 Pre SU Post SU 280 nm				●	●		
MCPXM20CNSBS2	Mobius® Chrom 20 Pre SU Post SU 280-300nm				●		●	
MCPXM20CNSBS3	Mobius® Chrom 20 Pre SU Post SU 254-280nm				●		●	
European Union								
MCPXM20EUSB01	Mobius® Chrom 20 Post SU 280 nm	Included				●		
MCPXM20EUSB02	Mobius® Chrom 20 Post SU 280-300nm						●	
MCPXM20EUSB03	Mobius® Chrom 20 Post SU 254-280nm							●
MCPXM20EUSBS1	Mobius® Chrom 20 Pre SU Post SU 280 nm				●	●		
MCPXM20EUSBS2	Mobius® Chrom 20 Pre SU Post SU 280-300nm				●		●	
MCPXM20EUSBS3	Mobius® Chrom 20 Pre SU Post SU 254-280nm				●		●	
Switzerland								
MCPXM20CHSB01	Mobius® Chrom 20 Post SU 280 nm	Included				●		
MCPXM20CHSB02	Mobius® Chrom 20 Post SU 280-300nm						●	
MCPXM20CHSB03	Mobius® Chrom 20 Post SU 254-280nm							●
MCPXM20CHSBS1	Mobius® Chrom 20 Pre SU Post SU 280 nm				●	●		
MCPXM20CHSBS2	Mobius® Chrom 20 Pre SU Post SU 280-300nm				●		●	
MCPXM20CHSBS3	Mobius® Chrom 20 Pre SU Post SU 254-280nm				●		●	
United Kingdom								
MCPXM20UKSB01	Mobius® Chrom 20 Post SU 280 nm					●		
MCPXM20UKSB02	Mobius® Chrom 20 Post SU 280-300nm						●	
MCPXM20UKSB03	Mobius® Chrom 20 Post SU 254-280nm							●
MCPXM20UKSBS1	Mobius® Chrom 20 Pre SU Post SU 280 nm				●	●		
MCPXM20UKSBS2	Mobius® Chrom 20 Pre SU Post SU 280-300nm				●		●	
MCPXM20UKSBS3	Mobius® Chrom 20 Pre SU Post SU 254-280nm				●		●	



Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97

Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

<https://millipore.nt-rt.ru> || <mailto:mer@nt-rt.ru>