

СЧЕТЧИК КЛЕТОК SCEPTER 2.0



Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97

Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

<https://millipore.nt-rt.ru> || <mailto:mer@nt-rt.ru>

What people are saying...

"At last, an alternative to lining up for the Coulter Counter®, and far easier than sweating over fragile hemocytometers."

AMY A. CAUDY is a Lewis-Sigler Fellow at Princeton University's Lewis-Sigler Institute for Integrative Genomics

The Scientist, Dec. 2010.
Top Ten Innovations of 2010.

"Cell counting is normally a very tedious process and usually only provides minimal information on the cell population. This instrument, which is only slightly larger than an automatic pipette, allows you to count cells in your tissue-culture hood, simplifies the procedure, and provides much useful data, such as the fraction of intact cells."

H. STEVEN WILEY is a lead biologist at the Environmental Molecular Sciences Laboratory at the Pacific Northwest National Laboratory

The Scientist, Dec. 2010.
Top Ten Innovations of 2010.



Scepter™ 2.0 – Precise, handheld cell counting

Scepter™ 2.0 is your portable cell counter. While other automated counters consume bench space and rely on object recognition software, manual focusing, and clumsy loading chambers, the Scepter™ cell counter provides true automation without the error that accompanies vision-based systems. With its microfabricated, precision-engineered sensor, the Scepter™ cell counter does all the work and delivers accurate and reliable cell counts in less than 30 seconds.

Scepter™ 2.0 marks the next generation in Scepter™ technology, highlighted by:

Compatibility with More Cell Types

The Scepter™ cell counter is the only one on the market to accurately count particles as small as 3 µm in diameter

Increased Cell Concentration Range

The new 40 µm sensor can count samples with concentrations as high as 1,500,000 cells/mL

Powerful Software for Complex, Effortless Cell Analysis

- Compare sample sets side by side using histogram overlay and multiparametric data table
- Create and save gating templates
- Generate reports, graphs and tables

The power of precision

Trust Scepter™ counting with your most valuable samples to get reproducible and reliable counts. The reliability of Scepter™ counting is particularly apparent with smaller cell types. Because the Scepter™ cell counter measures volume using the Coulter Principle, it can quantify cells based on size and will discriminate larger cells from smaller debris, unlike vision based techniques, which rely on object recognition software and cannot reliably detect small cells.

Scepter™ sensor technology

Compatible with 60 µm and 40 µm sensors, the Scepter™ 2.0 cell counter can meet even more of your cell- and particle-counting needs. Use the 60 µm sensor for particles between 6 and 36 µm. Use the 40 µm sensor for particles between 3 and 17 µm.

- Precise volumes are drawn into the Scepter™ sensor.
- As cells flow through the aperture in the sensor, resistance increases. This increase in resistance causes a subsequent increase in voltage.
- Voltage changes are recorded as spikes with each passing cell.
- Spikes of the same size are bucketed into a histogram and counted. This histogram gives you quantitative data on cell morphology that can be used to examine the quality and health of your cell culture.




Figure 1.

Particles are detected by Ohm's Law $V=IR$ (V =voltage, I =current, and R =resistance)

Cell Type	Measured size (µm)	40 µm sensor	60 µm sensor
2102 Ep	15-19		
454 beads			
A172	15		
A253	14-18		
A375	16		
A431	15-17		
A549			
Algae (various)	7-9		
B35	13-16		
B Cells	6-11		
C2C12	12		
C305	12-14		
C6	12-13		
CA46	10-12		
Caco-2	17		
CHO	14-17		
COS-1	12		
Cos-7	15		
D283	12		
Daudi	10-12		
DU-145	15-17		
Epithelia	14-15		
HCT-116	10		
HEK293	11-15		
HeLa	12-14		
HepG2	12		
HFF	18-20		
Hs27	14		
HT-1080	14-16		
HT-29	11		
HUH7- Hepatoma line			
Human ES Cells	9-12		
HUVEC	14-15		
IMR-32	12-14		
IMR-90	15		
Jurkat	13		
K562	22		
KB	14		
KG-1	10-13		
L6	14-16		
LNCaP	15-16		
Luminex® beads	5-6		
MCF7	15-17		
MDCK	13-15		
Meg-01	16-17		
MG-63	15-17		
Mouse ES Cell	5-13		
Mesenchymal Stem Cell	15-16		
MRC-5			
NCI-H146	10-13		
NIH 3T3	15		
NTERA2, clone D1	13		
OK	17-18		
PBMCs	7-12		
PC12	9-13		
Primary Astrocytes	7		
Primary Neuronal Cell			
Raji	12-15		
Ramos	11-12		
Rat Dorsal Root Ganglion Cells	7		
Rat Whole Blood	4.6		
Red Blood Cells	5-7		
Rat Neural Stem Cell	11-13		
RAW 264.7	12-15		
RBL	11-13		
RIN-mF5	13-14		
SF9	13		
SH-SY5Y	12		
Sk-Br-3	15-20		
SK-MEL-28	17-19		
SK-N-MC	14-15		
SK-N-SH	14-15		
Splenocytes	7-9		
SW-480	15		
SW-620	13-14		
T84	14-18		
T98G	17		
TF-1	13-14		
U251	16-20		
U20S	16-19		
U266	12		
U87-Human Glioblastoma cell line	12-14		
U937	11-13		
WI-38	12-15		
Y79	13-14		
Yeast- <i>Pichia Pastoris</i>	5		
Yeast- <i>S.cerevisiae</i>	6		

■ Recommended based on size ■ EMD Millipore Validated ■ Customer Validated

Table 1.

Cell types validated with the Scepter™ cell counter and the recommended Scepter™ sensor.

Scepter™ counting delivers precision

There is no need to subjectively determine cell counts, as required by vision-based counting methods. The Scepter™ cell counter detects every cell and displays the population as a histogram of cell size distributions. From the histogram, count all the cells or use the gating function to count a chosen subpopulation. By monitoring changes in your histogram, you can gain insight into the health and quality of your cell culture from one experiment to the next.





Figure 2. The average percent coefficient of variation (CV) for each counting method shown was calculated using cell concentration measurements at 50,000 cells/mL samples of 19 different cell lines. The Scepter™ cell counter is more precise than vision-based counting and hemocytometry, and approaches the precision of the Coulter Counter® standard (maroon bars). Error bars represent standard deviation.



	Format	Counting methods	Sample volume needed	Sample volume counted	Cells counted in a 100,000 cell/mL sample	Average % CV
Hemocytometer	Slide and microscope	Manual, vision-based	10 μL	.1 μL / square	10/square	41.8
Brand L	Benchtop	Automated vision-based system	10 μL	.4 μL	40	32.1
Scepter™ Cell Counter	Handheld	Impedance-based cell detection	100 μL	50 μL	5000	9.1

As easy as pipetting

Prepare the sample:

Start with a single-cell suspension, diluted to a total volume of 100 μL (recommended) in phosphate buffered saline (such as EmbryoMax® 1x DPBS) to 10,000-500,000 cells/mL (operating range for 60 μm sensor) in a 1.5 mL microcentrifuge tube.

Perform cell count:

- Turn on the Scepter™ cytometer by pressing the toggle on the back of the instrument and wait for on-screen instructions to appear.
- When prompted, attach a sensor to the end of the Scepter™ unit with the electrode sensing panel facing toward the front of the instrument, and you'll see detailed instructions for each step of the counting process.
- Pipette once to draw sample into the sensor. 50 μL of your cell suspension is drawn into the microfabricated, precision-engineered channel embedded in the sensor. The cell sensing zone detects each cell drawn into the sensor and thus cell concentration is calculated.
- The sensing zone also measures cell sizes and cell volumes with sub-micron and sub-picoliter resolution, enabling the Scepter™ cytometer to display a histogram distribution of cell size or cell volume.

Intuitive analysis software

From simple counts to complex volume measurements used to assess cell health parameters, Scepter™ Software Pro provides an intuitive, intelligent platform to perform high-level cell analysis based on the size measurements captured with the Scepter™ cell counter.

A VIEW OF SCEPTER™ SOFTWAREPRO



DATA:
data files from your
Scepter™ cell counter

CURRENT PLOT:
working plot and data file

GROUP STATS:
customizable statistics
from your selected data
files


MULTIFUNCTIONAL PLOT:
multiple data sets/
histogram overlays

REPORTS: export, print
selected graphs/files,
cut and paste

**ANALYSIS
TEMPLATES:** saved
gating parameters


Using the Scepter™ Software Pro on your computer, you can:

- Compare several samples and data sets side by side using histogram overlay and multiparametric tables
- Create and save gates to be used from one experiment to the next
- Create attractive graphical presentations and reports with your data



Ordering Information

Description	Qty	Catalogue No.
Scepter™ 2.0 Handheld Automated Cell Counter		
with 40 µm Scepter™ Sensors (50 Pack)	1	PHCC20040
with 60 µm Scepter™ Sensors (50 Pack)	1	PHCC20060
Includes:		
Scepter™ Cell Counter	1	
Downloadable Scepter™ Software	1	
O-Rings	2	
Scepter™ Test Beads	1	PHCCBEADS
Scepter™ USB Cable	1	PHCCCABLE
Scepter™ Sensors, 60 µm	50	PHCC60050
	500	PHCC60500
Scepter™ Sensors, 40 µm	50	PHCC40050
	500	PHCC40500
Universal Power Adapter	1	PHCCPOWER
Scepter™ O-Ring Kit, includes 2 O-rings and 1 filter cover	1	PHCCOCLIP



Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97

Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

<https://millipore.nt-rt.ru> || <mailto:mer@nt-rt.ru>